

Freely

Mobility Aids
Product Catalogue





Freely

Overview

Welcome to Freely and our range of mobility aids designed to enhance user independence and improve quality of life. The Freely range includes crutches, walking sticks, walkers and rollators to meet specific user needs and provide the utmost comfort and support.

At Freely we understand the significance of mobility assistance in maintaining an active and fulfilling lifestyle. Whether someone is recovering from an injury, living with a long-term condition, or simply seeking additional stability and confidence while on the move, our range of mobility aids has been carefully chosen to cater to their unique requirements.

From crutches that offer exceptional stability and weight distribution to walkers and rollators that combine durability with effortless manoeuvrability, each product in the range has been tried and tested to meet stringent safety standards and has been designed for the most important user requirements - ergonomics, lightweight materials, stability, anti-slip features, shock absorption, adjustability, portability and visual appeal.



Crutches

Elbow Crutch

An ergonomically designed elbow crutch to assist users with recovery, improve mobility and provide support for lower extremity injuries. The Elbow Crutch is sold in pairs and extensions are available (55cm in length).



Crutches

Permanent Crutch

The Permanent Crutch is a dependable solution for long-term walking support. The Permanent Crutch can be cut to size and is sold in pairs.



Colour Natural

Weight Limit

Standard 120kg
Heavy Duty 150kg

Height 120cm

**From cuff to ground*

Gutter Crutch

The Gutter Crutch is designed for individuals with arthritis, impaired grip, or partial weight-bearing conditions. This crutch features a padded forearm trough for even weight distribution providing stability for users. The adjustable hand grip ensures a comfortable and secure hold. This product is sold in single units.



Colour Natural

Adjustable Height

100-120cm

**From handle to ground*

Details

- Forearm padded with sponge and Velcro strap
- Adjustable handle from 10-20cm
- Height range from handle to ground 100-120,4cm

Crutches

Aluminium Axilla Crutches

The Aluminium Axilla Crutch is an underarm crutch designed for short term use as a result of a lower-limb injury and is available in five sizes. The Aluminium Axilla Crutch has soft rubber underarm padding and hand grips. Each size option is height adjustable and the handle position can also be adjusted for a more precise user fit. This product is sold in pairs.



Colour Natural

Weight Limit

XS-S 80kg
M-L 100kg
XL 114kg

Adjustable Height

XSmall 84-105cm
Small 94-115cm
Medium 115-135cm
Large 135-155cm
XLarge 145-165cm

**From underarm to ground*

Wooden Axilla Crutches

The Wooden Axilla Crutch is an underarm crutch designed for short term use as a result of a lower-limb injury and is available in three sizes. The Wooden Axilla Crutch has soft rubber underarm padding and hand grips. Each size option is height adjustable and the handle position can also be adjusted for a more precise user fit. This product is sold in pairs.



Colour Light Wood

Weight Limit

S 80kg
M 100kg
L 100kg

Adjustable Height

Small 96-117cm
Medium 110-130cm
Large 126-146cm

**From underarm to ground*

Walking Sticks

Classic Walking Stick

The perfect solution for the elderly, the physically challenged, or people who have trouble walking or standing. The Classic Walking Stick design is suitable for left and right handed users, and is sold in single units.



Weight Limit 110kg

Adjustable Height 70 - 90cm

**From handle to ground*



4 Colour Options

Quadruped Walking Stick

The perfect solution for the elderly, the physically challenged, or people who have trouble walking or standing. The design is suitable for left and right handed users, and is sold in single units.



Weight Limit 110kg

Adjustable Height 70 - 90cm

**From handle to ground*

Difference between the Classic and Quadruped Walking Stick

While the Classic Walking Stick relies on a single point of contact, the Quadruped Walking Stick offers four points of contact, providing increased stability and balance. The wider base of the Quadruped Walking Stick offers more support and helps distribute weight evenly, reducing the risk of falls or loss of balance..

Walking Sticks

Folding Walking Stick

The perfect solution for the elderly, those with mobility challenges, or difficulty walking. Suitable for both left and right-handed users, its foldable design adds convenience, making it easy to take anywhere, and is sold in single units.



Weight Limit 110kg

Folded Length 30cm

Adjustable Height 84-94cm

**From handle to ground*



2 Colour Options



S-Bend Walking Stick

The S-bend Walking Stick is a versatile aid designed to support walking and assist in maintaining good posture. This product is sold in single units.



Colour Natural

Weight Limit 100kg

Adjustable Height 72-94cm

Walking Pole

Designed for individuals transitioning from a walker, cane or crutch to a more independent gait, such as patients rehabilitating from injury or conditions like stroke. This product is sold in single units. Additional feet can be bought separately.



Colour Natural

Weight Limit 120kg

Standard Height 152cm

**Length can be cut to size*

Folding Walkers & Rollators

Folding Walker/Rollator Two Push Button

An easy-to-use folding walker designed as an ideal and dependable solution for those who require a supportive yet durable aid for their daily activities. This walker can be folded with the push of a quick double release button allowing for efficient storage and travel. The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.



Weight Limit 110kg

Adjustable Height

Adult 80 - 87cm

Youth 66 - 75cm

Paediatric 56 - 63cm

**From handle to ground*



Wheel Diameter 12,5cm



Folding Rollator

Wheels can be added to convert the walker into a rollator designed to be used primarily by people who are recuperating from leg or back injuries as well as hip and knee replacements. The two wheels at the front make for easier and more energy efficient propulsion whilst still providing a stable base of support.

Folding Walkers & Rollators

Reciprocal Folding Two Step Walker

An easy-to-use, adjustable Reciprocal Folding Walker with stand assist. The reciprocal walking frame walks with the user. There is no need to lift the walker forward with each step, rather this walker, when used correctly, steps with the user. This walker may also be used as a rigid frame.

This walker can be folded with the push of a quick single release button allowing for efficient storage and travel. The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.

Dimensions

Length 53cm, Width 47cm, Height 73 - 85cm



Weight Limit 110kg

Folding Walker One Push Button

An easy-to-use Folding Walker designed as an ideal and dependable solution for those who require a supportive yet durable aid for their daily activities. This walker can be folded with the push of a quick single release button allowing for efficient storage and travel.

The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.

Dimensions

Length 58cm, Width 49cm, Height 79 - 97cm



Weight Limit 110kg

Four Wheel Folding Rollator

A quality rollator designed to assist users with improved mobility. With an ergonomic design the lightweight frame is fitted with four large, sturdy castor wheels for efficient movement in every direction. Secure light touch hand brakes on both handles make stopping effortless, and both can be locked in place to ensure safety in a seated or standing position. The Rollator folds flat to a compact size for easy storage or transport.

Details

- Both handles are height adjustable
- Wide, comfortable seat and back rest (height 56cm)
- Standard seat (31.5 x 34.5cm), Wide seat (33.5 x 46.5cm)
- Detachable basket for storing belongings
- 18cm diameter wheels for both standard and wide

Weight Limit

Standard 100kg

Wide 136kg



Folding Walkers

Folding Walker (Youth)

An easy-to-use Folding Walker designed specifically for children who require a supportive yet durable aid for their daily activities. This walker can be folded with the push of a quick single release button allowing for efficient storage and travel.

The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.

Dimensions

Length 47cm, Width 58cm, Height 58 - 68cm



Folding Walker (Kiddies)

An easy-to-use Folding Walker designed specifically for children who require a supportive yet durable aid for their daily activities. This walker can be folded with the push of a quick single release button allowing for efficient storage and travel.

The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.

Dimensions

Length 39cm, Width 35cm, Height 49 - 54cm



Mobility Accessories

Toilet Seat Raiser

Used by patients especially after hip replacement surgery or other situations requiring an elevated seated position.

Details

- Fits all standard size toilets
- Made from white plastic (LLDPE)
- Features stainless steel feet with rubber tips
- Non-adjustable



Standard Grey Cuffs

Designed for the standard adult elbow crutch these replacement cuffs are sold in pairs.

Colour Grey



Standard Crutch Feet

Standard crutch feet made from grey PVC, a thick base and galvanised washer inside. The standard Crutch Feet are available in 19, 21 and 27mm diameter size options.

Colour Grey

Details

- Should be replaced regularly
- Easy to attach



Partner with Us

At Freely Mobility, we take pride in being a trusted supplier of mobility aids. We understand the importance of providing you with quality products at competitive prices that will positively impact your business. Our goal is to help you build a reputation for delivering reliable mobility aid products to your own customers, enhancing customer satisfaction, and overall competitiveness in the market.

Here are some key advantages of partnering with us:

Strong Partnerships

We prioritise building long-term partnerships fostering open communication, trust, and mutual growth.

Expertise and Experience

We have years of experience sourcing, manufacturing and supplying quality products throughout Southern Africa.

Wide Product Range

Our products cater to various end user needs, providing you with flexibility in your offering and we're constantly developing or sourcing new, innovative products.

High-Quality Products

We provide reliable and durable products ensuring consistent product quality for your customers.

Freely

 www.freelymobility.co.za

 info@freelymobility.co.za